

# Beyes®

The Next Standard

## SHOW PROMOTIONS

Only valid during the show Through Oct 6

ODA Booth #336 : [www.beyes.ca/en/events](http://www.beyes.ca/en/events)



### AirLight® M800 Plus

LED High Speed Handpiece

- 23W of Torque
- Ceramic Ball Bearings
- Handles up to 60 PSI of air pressure
- Titanium Body

| Backend | Code    |
|---------|---------|
| 4 Hole  | HP3008P |
| Beyes®  | HP6018P |
| W&H®    | HP3022P |
| Star®   | HP3054P |
| NSK®    | HP3043P |
| KAVO®   | HP3036P |

Reg. \$498/ea

**Sale: 3 For \$1000 Only**



### Slimax-C Plus®

LED Curing Light System

- Built-In Digital Radiometer
- LED Display Panel
- Cordless, Portable Design
- Cures 850 Times Per Charge
- Up to 2000 mW/cm ± 10%
- Rechargeable Samsung Lithium-Ion Battery
- Solid Aluminum Body
- 1 Year Warranty

| Color  | Code       |
|--------|------------|
| Silver | SM1003P-S  |
| Blue   | SM1003P-B  |
| Pink   | SM1003P-PI |

Reg. \$598/ea

**Buy 1 Get 1 FREE**

**Net = \$299/ea**



### Maxso® E600

Electric Micromotor System

This newest brushless electric micromotor system is designed with ultimate user-friendliness in mind. Its solid titanium construction ensures perfect balance - it is also one of the most lightweight systems of its kind available. Integrate the Maxso E600 into your practice with ease and start experiencing the advantages of electric today.

- **Powerful:** 3.0 Ncm of torque
- **Steady Performance:** Speed of rotations is controlled between 2,000 and 40,000min-1
- **Smooth Operation:** Torque ripple and heat reduction technology eliminates vibrations
- **Convenient & Durable:** Autoclavable
- **Easy Installation:** Simply mount and connect



Code: MT6001 Reg. \$1,598

#### Package 1:

**\$1,999 (Savings of: \$543)**

Includes: One E600, One X99L Attachment



#### Package 2:

**\$2,599 (Savings of: \$887)**

Includes: One E600, Two X99L Attachments



#### Attachments:

### Maxso® X99L

Contra Angle

- 1:5 Ratio
- Anti-Heat System
- 4 Points Spray
- Ceramic Ball Bearing
- Fiber-optic
- 12 Months Warranty

Code: SL3025

Reg. \$944/ea

**Buy 4 Get 1 FREE**

**Net = \$755.20/ea**



### Maxso® M29-MH01

Contra Angle

- 1:1 Ratio
- Internal Spray
- Fiber-Optic
- Push Button
- For CA Burs/Latch Type
- 12 Months Warranty

Code: SL3068

Reg. \$355/ea

**Buy 4 Get 1 FREE**

**Net = \$284/ea**



### Maxso® S20A-IS

Straight Nose Cone

- 1:1 Ratio
- Internal Spray
- Fiber-Optic
- For HP bur
- 12 Months Warranty

Code: ST2009

Reg. \$376/ea

**Buy 4 Get 1 FREE**

**Net = \$300.80/ea**



## Airlight® M800X

LED High Speed Handpiece

- Instant-Stop Capability
- 27W of Torque
- X-Ball Bearings
- Extra Quiet
- Micro-Tex Slip-Proof Coating
- 25 Month Warranty

| Backend | Code    |
|---------|---------|
| Beyes®  | HP6018X |
| Star®   | HP3054X |
| NSK®    | HP3043X |
| KAVO®   | HP3036X |

Reg. \$1,198/ea  
**Buy 3 Get 1 FREE**  
Net = \$898.50/ea



## Airlight® M800 Surgical 45

LED High Speed Handpiece

- 45 Degree Angled Head
- Rear Air Exhaust
- Direct LED
- Stream Water
- 22W of Cutting Power

| Backend | Code   |
|---------|--------|
| 4 Hole  | HP3048 |
| Beyes®  | HP6038 |
| NSK®    | HP3044 |
| KAVO®   | HP3037 |

Reg. \$560 **Sale: \$498/ea**



## Maxso® X200 Smart

Fiber Optic High Speed Handpiece

- Triple Water Spray
- Push Button Chuck
- 23W of Torque
- Ceramic Ball Bearings
- Swivel

| Backend | Code   |
|---------|--------|
| KAVO®   | HP2060 |
| NSK®    | HP2066 |

Reg. \$438/ea  
**Buy 3 Get 1 FREE**  
Net = \$328.50/ea



## Maxso® M200E

Non Optic High Speed Handpiece

- Single Spray
- Push Button
- Stainless Steel Ball Bearings
- Low Turbine Noise
- No Swivel

| Backend | Code   |
|---------|--------|
| 4 Hole  | HP2105 |

Reg. \$106/ea  
**Buy 5 Get 1 FREE**  
Net = \$88.33/ea



## CareMaster®

Automatic Handpiece Maintenance System

Maintain Your Most Significant Investments in the Office!  
Clean 3 Handpieces Per Cycle!

- Maximize the lifespan of handpieces
- Economical, efficient and time-saving
- Filtered air used to rotate the handpieces
- Detachable door for easy cleaning
- Provides care for all brands of handpieces



Code: HPL2010  
2 Turbines & 1 E-Type Attachment,  
**\$1,154**

Code: HPL2011  
1 Turbines & 2 E-Type Attachments,  
**\$1,258**

\* Requires Airline Connection  
\* Coupler Sold Separately

Buy Any Model of CareMaster  
**Get 1 AirLight® M800 Plus FREE**  
(\$498 Value)

## EZ-Swivel®

Couplers



Beyes® PD  
4 Hole  
Code: HP6008  
**\$136**



Kavo® Type  
4 Hole, Non-optic  
Code: HP2063  
**\$146**



Kavo® Type  
6 Hole, Fiber-optic  
Code: HP2064  
**\$156**

## Maxso® X92

Implant Handpiece

- 20:1 Ratio
- Non-Optic
- Push Button
- CA Burs
- Fits into major implant motors

Code: SL3019

Reg. \$799 **Sale: \$499**



**Maxso®**

Low Speed Package

Includes:  
Contra Angle-Push Button (SL2026)  
Straight Nose Cone (ST2010)  
Air Motor-NS/M4 Backend (MT2010)



~~Reg. \$502/pk~~

**Sale: \$399/pk**

**EP12-NS**

Prophy Handpiece

- 4:1 Ratio
- For disposable prophy angles



Code: SL1002

~~Reg. \$198/ea~~

**Buy 4 Get 1 FREE**

**Net = \$158.40/ea**

**Maxso®**

Low Speed Attachments

- 1:1 Ratio
- E-Type

**Buy 3 Get 1 FREE**

**Snap-On**

Code: SL2033

~~Reg. \$104/ea~~

**Net = \$78/ea**



**Screw-On**

Code: SL2034

~~Reg. \$104/ea~~

**Net = \$78/ea**



**Latch Type, CA Burs/Latch**

Code: SL2028

~~Reg. \$114/ea~~

**Net = \$85.50/ea**



**Push Button, CA Burs/Latch**

Code: SL2026

~~Reg. \$146/ea~~

**Net = \$109.50/ea**



**Straight Nose Cone, HP Burs**

Code: ST2010

~~Reg. \$94/ea~~

**Net = \$70.50/ea**



\* All items (purchased and free) under the 3+1 promotion must be of the same model/code. No mixing and matching.

**Maxso®**

Head

**Buy 4 Get 1 FREE**

**Snap-On**

Code: SL2015

~~Reg. \$46/ea~~

**Net = \$36.80/ea**



**Screw-On**

Code: SL2016

~~Reg. \$46/ea~~

**Net = \$36.80/ea**



**Latch Type**

Code: SL2010

~~Reg. \$52/ea~~

**Net = \$41.60/ea**



**Push Button**

Code: SL2008

~~Reg. \$90/ea~~

**Net = \$72/ea**



**Maxso®**

Air Motor

E-Type, No Spray, 20k, 4 holes

Code: MT2010

~~Reg. \$262/ea~~

**Buy 4 Get 1 FREE**

**Net = \$209.60/ea**



**Endopilot®**

Cordless Endodontic Handpiece

- Operates continuously for more than 2 hours
- Memory for up to 5 programs.
- Charging interface is designed as MICRO USB.

Code: AL2030

~~Reg: \$999~~

**Sale: \$882**



**Apexipilot® All-In-One**

Endodontic Handpiece, Built-in Apex Locator

- Built-in apex locator
- Quick start mode
- Bluetooth

Code: AL2020

~~Reg: \$1,699~~

**Sale: \$1,508**



**Apexipilot®**

Apex Locator

- Automatic calibration ensuring accurate measurement
- Bluetooth wireless transmission
- Autoclavable file clip, lip hook and touch probe
- Built-in non-removable battery

Code: AL2040

~~Reg: \$649~~

**Sale: \$499**



**Accurator®**

Apex Locator

Measure the length of canals in all conditions including: dry, wet and where blood is present with industry-leading accuracy.

- Automatic Calibration
- Autoclavable Accessories

Code: AL2010

~~Reg: \$549~~

**Sale: \$349**



## Comfort-Sonic®

### LED Piezo Scaler

- 3 Modes: Scaling, Perio and Endo
- Maximizes Visibility
- Minimizes Eye Strain
- Detachable Handpiece
- Autoclavable Up to 135 °C

- Best illumination - 6 LEDs through an optical ring
- Comfort sonic engine - stable power output and performance
- Extensive tip selection for multiple applications
- Built-in solution container (for P9)
- Integrated water pump (for P9)



#### P6 Plus

Code: UL2104  
6 tips included

Reg. 2 for \$999

**Sale: 2 For \$800 Only**

#### P9

Code: UL2103  
8 tips included, 2 water bottles

Reg. \$1,398.95

**Sale: \$699**

## Comfort-Sonic®

### Piezo Scaler Handpiece & Accessories

#### Non-optic:

|                      |             |
|----------------------|-------------|
| <b>Compatibility</b> | <b>Code</b> |
| Beyes & EMS          | UL2020      |
| Satelec®             | UL2021      |

Reg. \$180.84/ea

**Buy 3 Get 1 FREE**

**Net = \$135.63/ea**



#### Fiber-optic:

|                      |             |
|----------------------|-------------|
| <b>Compatibility</b> | <b>Code</b> |
| Beyes & EMS          | UL2022      |
| Satelec®             | UL2023      |

Reg. \$212.94/ea

**Buy 3 Get 1 FREE**

**Net = \$159.71/ea**



#### Tip Wrench:

WST

Code: UL2420

~~REG: \$17.12~~

**Sale: \$16<sup>99</sup>**



## Canaview®

### Intraoral Camera

- True Auto-Focus
- Innovative Anti-Fog Design
- Responsive Capture Button
- 5 MP High Resolution Images
- Intuitive Camera Lighting System
- No Manual Adjustments Required
- Direct USB Connection to Computer
- Super Light-Weight
- Complimentary Image Management Software

Code: EV2060

Reg. \$838/ea

**Buy 3 Get 1 FREE**

**Net = \$628.50/ea**



#### Scaling Tips:

Used to remove calculus and plaque on supragingival, neck of teeth and hard calculus on the teeth point.



|              |             |
|--------------|-------------|
| <b>Model</b> | <b>Code</b> |
| E-G1         | UL2300      |
| S-G1         | UL2400      |

Used to remove all supragingival calculus and plaque.



Equivalent to EMS\* Type A

|              |             |
|--------------|-------------|
| <b>Model</b> | <b>Code</b> |
| E-G4         | UL2302      |
| S-G4         | UL2402      |

#### Perio Tips:

Used to remove subgingival calculus.



|              |             |
|--------------|-------------|
| <b>Model</b> | <b>Code</b> |
| E-P1         | UL2305      |
| S-P1         | UL2405      |

Used to remove subgingival and deep perio calculus.



Equivalent to EMS\* Type PS

|              |             |
|--------------|-------------|
| <b>Model</b> | <b>Code</b> |
| E-P3         | UL2306      |
| S-P3         | UL2406      |

Reg. \$29.95/ea

**Promo 1:**

**Buy 5 @ \$25.96 each**

**Promo 2:**

**Buy 10 @ \$21.96 each**

Beyes Dental Canada Inc.

23-595 Middlefield Rd., Toronto, ON, M1V 3S2, Canada | T: 1-855.603.1888 | F: 1-855.720.1228 | E: order@beyes.ca

\*Beyes Dental Canada Inc. is not responsible for any typographical errors. Specifications are subject to change without prior notice. All brands are holders of their respective trademarks. File # ENF010R1

Hughes Diversified Dental LLC 800-773-0800

**Beyes®**  
The Next Standard